Chartis®

Pulmonary Assessment System Instructions for Use

STERILE EO

1.0 Product Description

The CHARTIS® Pulmonary Assessment System consists of two components and is designed to be used by experienced bronchoscopists during a diagnostic bronchoscopy.

1. Chartis Catheter and Connector Set
2. Chartis Console, (see User Manual for Instructions for Use of the Console)

The Chartis Catheter (Figure 1) is a single use, sterile, disposable device designed to be inserted into a 2.8mm working channel of a standard video bronchoscope during a diagnostic bronchoscopy procedure. After the target lung segment is accessed by the bronchoscope, the distal tip of the Chartis Catheter can be pushed out from the distal end of the bronchoscope directly into the target airway. Inflation of the compliant balloon component on the distal tip of the Chartis Catheter causes the airway to become sealed. Air can then flow out from the target compartment into the environment only through the Chartis Catheter central lumen. Airway flow and pressure exiting the Chartis Catheter during spontaneous respiration is displayed in a graphic format on the Chartis Console. Functional assessment of air pressures and flows takes less than 5 minutes in isolated lung compartments. After assessment the catheter is removed from the airway.

The Chartis Catheter is a hollow catheter shaft made of medical-grade PEBAX. The Chartis Catheter is 169 cm in length with a 2.7 mm outer diameter. The distal end of the catheter employs an inflatable balloon made from a compliant elastomeric material. The catheter has two lumens: the main lumen allows air flow from the isolated lung compartment to the Console, and the second smaller lumen allows balloon inflation. The proximal end of the catheter has two ports, a balloon inflation port to inflate the distal balloon, and an obturator port with a luer fitting for connection to the Chartis Console via the Connector Set. The design of the balloon is similar to other compliant bronchial balloons. A custom obturator inserted into the main lumen of the catheter is used for both pushability and to clear mucus that may obstruct the lumen during catheter manipulation in the airways. The only portion of the device that comes into bodily contact is the distal tip of the catheter, which is made from medical grade PEBAX (shaft) and a compliant elastomeric material (balloon). The Chartis Catheter is designed to capture information at the lobar and segmental level. No medicinal substances are eluted from the catheter.



Figure 1: Pulmonx Chartis Catheter and Distal Balloon Tip

Components:

- Chartis Catheter
- Connector Set

Specifications:

Catheter OD (bronchoscope interfacing)	2.7 mm
Working Length (CHR-CA-12.0)	72 cm
Working Length (CHR-CA-12.0-XL)	76 cm
Total Length	169 cm
Airway Diameter Range	5mm-12 mm
Maximum air volume	3 ml (cc)

2.0 Indications for Use

The Chartis System is indicated for use by bronchoscopists during a diagnostic bronchoscopy in adult patients with Chronic Obstructive Pulmonary Disease (COPD) and emphysema in a bronchoscopy suite. The system, composed of the Chartis Catheter and Chartis Console, is designed to measure pressure and flow in order to calculate resistance to airflow and quantify collateral ventilation in isolated lung compartments. The Chartis Catheter is used through the working channel of a bronchoscope and connects to the Chartis Console. The Chartis Console is a reusable piece of capital equipment that displays the patient information.

3.0 Contraindications

The Chartis System is contraindicated in the presence of active infection or major bleeding diathesis. There are no known interfering substances.

4.0 Warnings and Precautions

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to patient injury.
- The device should only be used by physicians trained in interventional bronchoscopic techniques.
- The Chartis Catheter is designed for use in bronchoscopes with a minimum working channel of 2.8 mm.
- The Chartis Catheter is designed to be used only in conjunction with the Chartis Console.
- Use of the Chartis Catheter in the main stem bronchus is not recommended.
- For Standard Mode, the patient must have spontaneous respiration during the air flow and pressure measurements. For Ventilator Mode, refer to the Chartis Console User Manual for further information.
- Over tightening of luer fittings may cause damage to the device.
- Under tightening of luer fittings may cause an inaccurate assessment.
- Balloon inflation in an improper location may lead to patient injury.
- Do not exceed a maximum inflation volume of 3 mL (cc).
- Over inflation of balloon may lead to balloon burst.
- Repositioning of catheter while balloon is inflated may result in deformation of balloon and inaccurate assessment.
- If the balloon bursts during use, remove the Chartis Catheter immediately and discard.
- For **CHR-CA-12.0**: Avoid disengaging the catheter strain relief from the bronchoscope valve while the obturator is removed from the catheter to prevent kinking.
- For **CHR-CA-12.0-XL**: To prevent kinking, avoid excessive bending of the catheter at the bronchoscope valve while the obturator is removed from the catheter.
- The Chartis Catheter is intended for single-patient use only. Do not re-sterilize. No assurance of sterility can be made if devices are re-used. Do not use the device if the sterilization barrier has been damaged or if the device is dropped after removal from sterile packaging. Re-use of device or use of non-sterile device may result in patient infection
- Do not use if the package or contents appear to be damaged. Use aseptic technique when handling.
- Use of the Chartis Catheter when the patient is moved into the lateral decubitus position may result in difficulty in the initial placement of the catheter.

5.0 Potential Complications

Potential complications include, but are not limited to, the following:

- Allergy to catheter materials
- Epistaxis
- Aphonia
- Fever
- Bleeding
- Headache
- Bronchospasm
- Hemoptysis
- Chest pain
- Hypotension
- COPD exacerbation
- Hypercapnea
- Cough
- Hypoxemia
- Disorientation / anxiety
- Impaired lung function
- Embolism
- Increased mucus secretions

- Infection
- Pneumothorax
- Irritation at balloon site
- Sepsis
- Iatrogenic injuries
- Self Limiting Arrhythmia
- Nausea / vomiting
- Shortness of breath
- Pain
- Sore throat
- Perforation
- Tissue irritation
- Pneumonia
- Vocal cord injury

6.0 Required Accessories

- 3cc syringe, luer-lock, sterile, for balloon inflation
- 10cc syringe, luer-lock, sterile, for catheter flushing.
- One-way stopcock, sterile (female luer lock to male swivel connector)
- Three-way stopcock sterile (two female luer locks and one male swivel connector)

7.0 Procedure

7.1 Preparation

1. Inspect sterile packaging for any signs of damage. Do not use if the packaging appears to be damaged.
2. Open the Chartis Catheter and Console Connector Set onto a sterile field for preparation.
3. Inspect the devices for any damage.
4. Retract the plunger on the inflation syringe and attach the syringe and the one-way stopcock to the balloon inflation port (labeled "BALLOON") on the proximal portion of the Chartis Catheter
5. Test the balloon by injecting approximately 1 mL (cc) of air. Inspect the balloon to ensure that it holds air volume. If the balloon holds air volume it is ready for use. Retract the plunger of the syringe to remove the injected air.
6. Ensure the Chartis Console is set-up and prepared in accordance with the Console User's Manual.
7. Inspect the bronchoscope for any damage (e.g. sharp edges or cracks) to prevent damage to the Chartis Catheter.
8. Ensure bronchoscope is clean and free of particulate to prevent damage to the Chartis Catheter.

7.2 Usage

1. Inspect the airway during a standard diagnostic bronchoscopy procedure. Use suction to remove any excessive mucus.
NOTE: The patient must have spontaneous respiration during the air flow and pressure measurements.
2. Attach the Connector Set to the Chartis Console.
3. Lubricate the distal tip of the Catheter with a viscous sterile lubricant for ease in passing the Catheter through the working channel of the bronchoscope.
4. Attach one way stopcock to balloon inflation port of the catheter.
5. Depress plunger of sterile 3cc syringe and attach to opposite end of stopcock. Ensure stopcock is in open position.
6. Apply a vacuum to the balloon lumen by retracting the plunger of the inflation syringe and closing the stopcock prior to inserting the catheter into the bronchoscope.

7. Carefully insert the Chartis Catheter with the obturator into the working channel of the bronchoscope.
8. Position the bronchoscope into the target airway (e.g., RUL B1).
9. For **CHR-CA-12.0**: Advance the Chartis catheter approximately 2 cm beyond the distal tip of the bronchoscope until the catheter's distal balloon tip is visualized and the strain relief is engaged in the bronchoscope valve.

NOTE: Use caution to avoid kinking the catheter when the strain relief marker is visible outside the bronchoscope valve.

Figure 2: For CHR-CA-12.0 the Visible Strain Relief Marker indicates catheter not fully inserted

10. For **CHR-CA-12.0-XL**: Advance the Chartis catheter approximately 2 cm beyond the distal tip of the bronchoscope until the catheter's distal balloon tip is visualized.

NOTE: Use caution to avoid kinking the catheter at the point where it is inserted into the bronchoscope valve.
11. Position the distal tip of catheter into the target airway. Small positioning adjustments with the bronchoscope and/or the Chartis Catheter should be made to ensure the balloon is positioned in the target airway location. The black proximal balloon marker identifies the proximal edge of the balloon.
12. Remove the obturator from the catheter.
13. Connect the three-way stopcock to the Connector Set and connect the other end to the main lumen port on the catheter handle.
14. Flush the catheter.
 - Retract the plunger on a 10cc syringe and connect it to the third port on the three-way stopcock.
 - Position the three-way stopcock in such that the syringe port is in the "OFF" position.
 - Depress plunger of the 10cc syringe to build up pressure.
 - Rotate the stopcock valve to allow flow from the syringe to the catheter. Do not direct flow to the console.
 - Detach the 10cc syringe from the stopcock and repeat the previous steps above to flush the catheter a second time.

- Once the flushing process is complete turn the stopcock valve such that the syringe port is in the "OFF" position to ensure a free flow from the catheter to the console.

NOTE: Do not apply vacuum to the catheter or console during the flushing process.

15. Using the inflation syringe, inject up to 3 mL (cc) of air as needed for balloon expansion to create a seal in the airway under constant visualization, and close the stopcock to secure the inflated balloon. Do not over-inflate balloon in smaller airways.

CAUTION: To prevent balloon rupture, do not inject more than 3 ml (cc) of air.

16. Confirm via visual inspection that the target airway is occluded. Air bubbles at the balloon should not be visible during respiration.
17. The system is now ready for data capture. See the Operators Manual for the Chartis Console for complete instructions for use.
18. Hold the bronchoscope/Chartis catheter in place while the Chartis Console takes measurements.

NOTE: If the Chartis Catheter becomes blocked or clogged during use, insert the obturator to clear the lumen or small amounts of air may be injected through the inner lumen of the Catheter as described in the "Flush the catheter" section above.

19. When a measurement is complete (normally less than 5 minutes at each isolated lung segment), open the stopcock and retract the plunger on the inflation syringe to deflate the balloon on the Chartis Catheter.
20. Re-insert the obturator.
21. Retract the Chartis Catheter into the distal tip of the bronchoscope.
22. Reposition the bronchoscope near the next segment or lobe and repeat steps above to assess other airways as required.
23. When all measurements are complete, deflate the balloon, insert obturator and remove the Chartis Catheter by retracting it through the working channel of the bronchoscope.
24. Dispose of the Chartis Catheter per institutional standard for biohazardous materials.
25. Any serious incidents or serious adverse events that occur due to the use of the Chartis Catheter should be reported to Pulmonox (the manufacturer) and any relevant regulatory bodies (e.g., the European Competent Authority in the relevant EU Member State).

8.0 How Supplied

- The Chartis Catheter is supplied sterile in a peel-open package.
- The device is intended for single-patient use only. Do not re-sterilize the device.
- Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken.




9.0 Storage and Handling

- Storage: Store at room temperature in a clean and dry environment. Keep dry.
- Protect Catheter from UV light exposure.
- Handle with care
- Avoid extremes of temperature and humidity

10.0 Sterile

Sterilized with ethylene oxide. Do not use if the package is damaged or open.

Graphic Symbols Contained in Device Labeling

	Batch Code		Do Not Re-sterilize
	Catalog Number		Do Not Reuse
	Contents		Manufacturer
	Non-Pyrogenic		Date of Manufacture
	Sterilized Using Ethylene Oxide		Consult Instructions for Use
	Use By		Keep Dry
	Medical Device		Unique Device Identifier
	Single Sterile Barrier System		Importer
	Caution, consult accompanying documents		Authorized Representative in the European Community/EU
	Authorized Representative in Switzerland/CH		Authorized Representative in the United Kingdom UK
	Temperature limit		
		Maximum Inflation Volume	
	Do Not Use if Package is Damaged or Opened		
	Minimum Brochoscope Working Channel Diameter		

This product and/or its use are covered by one or more of the following United States patents: 7,883,471; 8,342,182; 8,523,782; 8,454,527. Other U.S. patents pending.

This product and/or its use are covered by one or more of the following international patents: EP07841495; DE6207028454.6; GB2056715; FR2056715; EP02768834; AT1435833; BE1435833; CH1435833; DE1435833; FR1435833; GB1435833; IE1435833; IT1435833; NL1435833; TR1435833; DE2056715; GB2056715; FR2056715; JP2008528105. Other international patents pending.

 **Pulmonx Corporation**
700 Chesapeake Drive
Redwood City, CA 94063
United States of America
Tel: +1 (650) 364-0400

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 **MDSS-UK RP Ltd.**
6 Wilmslow Road, Rusholme
Manchester M14 5TP
United Kingdom

 **Pulmonx International Sàrl**
Rue de la Treille 4, 2000
Neuchâtel, Switzerland
Tel: +41 32 475 2070

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Pulmonx GmbH