Chartis® Precision XL Catheter Pulmonary Assessment System

Instructions for Use

1.0 PRODUCT DESCRIPTION

The CHARTIS® Pulmonary Assessment System consists of two components and is designed to be used by experienced bronchoscopists during a diagnostic bronchoscopy.

1. Chartis Precision XL Catheter, Connector Set, 3 mL Syringe, and 10 mL Syringe
2. Chartis Console, (see User Manual for Instructions for Use of the Console)

The Chartis Precision XL Catheter (Figure 1) is a single use, sterile, disposable device designed to be inserted into a 2.8mm working channel of a standard video bronchoscope during a diagnostic bronchoscopy procedure. After the target lung segment is accessed by the bronchoscope, the distal tip of the Chartis Precision XL Catheter can be pushed out from the distal end of the bronchoscope directly into the target airway. Inflation of the compliant balloon component on the distal tip of the Chartis Precision XL Catheter causes the airway to become sealed. Air can then flow out from the target compartment into the environment only through the Chartis Precision XL Catheter central lumen. Airway flow and pressure exiting the Chartis Precision XL Catheter during spontaneous respiration is displayed in a graphic format on the Chartis Console. Functional assessment of air pressures and flows takes less than 5 minutes in isolated lung compartments. After assessment the catheter is removed from the airway.

The Chartis Precision XL Catheter is a hollow catheter shaft made of medical-grade PEBAX. The Chartis Precision XL Catheter is 145 cm in length with an outer diameter less than 2.8 mm. The distal end of the catheter employs an inflatable balloon made from a compliant elastomeric material. The catheter has two lumens: the main lumen allows air flow from the isolated lung compartment to the Console, and the second smaller lumen allows balloon inflation. The proximal end of the catheter has two ports, a balloon inflation port to inflate the distal balloon, and an obturator port with a luer fitting for connection to the Chartis Console via the Connector Set. The design of the balloon is similar to other compliant bronchial balloons. A custom obturator inserted into the main lumen of the catheter is used for both pushability and to clear mucus that may obstruct the lumen during catheter manipulation in the airways. The only portion of the device that comes into bodily contact is the distal tip of the catheter, which is made from medical grade PEBAX (shaft) and a compliant elastomeric material (balloon). The Chartis Precision XL Catheter is designed to capture information at the lobar and segmental level. No medicinal substances are eluted from the catheter.



Figure 1: Pulmonx Chartis Precision XL Catheter and Distal Balloon Tip

Components

- Chartis Precision XL Catheter
- Connector Set
- 3 mL Syringe
- 10 mL Syringe

Specifications

Catheter OD (bronchoscope interfacing)	<2.8 mm
Working length	72 cm
Total length	145 cm
Balloon Diameter Size Range (to maximum inflation volume)	5 mm – 13 mm
Maximum air volume	4 mL (cc)

2.0 INDICATIONS FOR USE

The Chartis System is indicated for use by bronchoscopists during a bronchoscopy in adult patients with Chronic Obstructive Pulmonary Disease (COPD) and emphysema in a bronchoscopy suite. The system, composed of the Chartis Catheter and Chartis Console, is designed to measure pressure and flow in order to calculate resistance to airflow and quantify collateral ventilation in isolated lung compartments. The Chartis Catheter is used through the working channel of a bronchoscope and connects to the Chartis Console. The Chartis Console is a reusable piece of capital equipment that displays the patient information.

3.0 CONTRAINDICATIONS

The Chartis System is contraindicated in the presence of active infection or major bleeding diathesis. There are no known interfering substances.

4.0 WARNINGS AND PRECAUTIONS

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to patient injury.
- The device should only be used by physicians trained in interventional bronchoscopic techniques.
- The Chartis Precision XL Catheter is designed for use in bronchoscopes with a minimum working channel of 2.8 mm.
- The Chartis Precision XL Catheter is designed to be used only in conjunction with the Chartis Console.
- Use of the Chartis Precision XL Catheter in the main stem bronchus is not recommended.
- For Standard Mode, the patient must have spontaneous respiration during the air flow and pressure measurements. For Ventilator Mode, refer to the Chartis Console User Manual for further information.
- Over tightening of luer fittings may cause damage to the device.
- Under tightening of luer fittings may cause an inaccurate assessment.
- Balloon inflation in an improper location may lead to patient injury.
- Do not exceed a maximum inflation volume of 4 mL (cc).
- Over inflation of balloon may lead to balloon burst.
- Repositioning of catheter while balloon is inflated may result in deformation of balloon and inaccurate assessment.
- If the balloon bursts during use, remove the Chartis Precision XL Catheter immediately and discard.
- To prevent kinking, avoid excessive bending of the catheter at the bronchoscope valve while the obturator is removed from the catheter.
- Upon removal from the packaging and sterile pouch, protect the balloon material from excessive exposure to UV radiation.
- For single-patient use only.
- Do not re-sterilize. No assurance of sterility can be made if devices are re-used. Do not use the device if the sterilization barrier has been damaged or if the device is dropped after removal from sterile packaging. Re-use of device or use of a non-sterile device may result in patient infection.
- Do not use if the package or contents appear to be damaged. Use aseptic technique when handling.
- Use of the Chartis Precision XL Catheter when the patient is moved into the lateral decubitus position may result in difficulty in the initial placement of the catheter.
- The Chartis System safety and performance will be made available and will be kept up to date on Pulmonx's website: www.pulmonx.com/safety-performance/

5.0 COMPLICATIONS

Potential complications include, but are not limited to, the following:

- | | |
|---------------------------------|----------------------------|
| • Allergy to catheter materials | • Cough |
| • Aphonia | • Disorientation / anxiety |
| • Bleeding | • Embolism |
| • Bronchospasm | • Epistaxis |
| • Chest pain | • Fever |
| • COPD exacerbation | • Headache |

- | | |
|------------------------------|----------------------------|
| • Hemoptysis | • Pain |
| • Hypotension | • Perforation |
| • Hypercapnea | • Pneumonia |
| • Hypoxemia | • Pneumothorax |
| • Impaired lung function | • Sepsis |
| • Increased mucus secretions | • Self Limiting Arrhythmia |
| • Infection | • Shortness of breath |
| • Irritation at balloon site | • Sore throat |
| • Iatrogenic injuries | • Tissue irritation |
| • Nausea / vomiting | • Vocal cord injury |

6.0 REQUIRED ACCESSORIES

- 3 mL syringe, luer-lock, sterile, for balloon inflation
- 10 mL syringe, luer-lock, sterile, for catheter flushing

7.0 PROCEDURE

7.1 Preparation

1. Inspect sterile packaging for any signs of damage. Do not use if packaging appears damaged.
2. Open the Chartis Precision XL Catheter and Console Connector Set onto a sterile field for preparation.
3. Inspect the devices for any damage.
4. Retract the plunger on the 3 mL inflation syringe.
5. Attach the 3 mL syringe to the balloon inflation port of the one-way stopcock on the handle of the Chartis Precision XL Catheter. See Figure 2 below.


Figure 2: Balloon Inflation Port Connection

6. Test the balloon by injecting approximately 1 mL (cc) of air. Inspect the balloon to ensure that it holds air volume. If the balloon holds air volume, it is ready for use. Retract the plunger of the syringe to remove the injected air.
7. Ensure the Chartis Console is set-up and prepared in accordance with the Console User Manual.
8. Inspect the bronchoscope for any damage (e.g., sharp edges or cracks) to prevent damage to the Chartis Precision XL Catheter.
9. Ensure bronchoscope is clean and free of particulate to prevent damage to the Chartis Precision XL Catheter.

7.2 Usage

1. Inspect the airway during a standard diagnostic bronchoscopy procedure. Use suction to remove any excessive mucus.

NOTE: For Standard Mode, the patient must have spontaneous respiration during the air flow and pressure measurements. For Ventilator Mode, refer to the Chartis Console User Manual for further information.

2. Attach the Connector Set to the Chartis Console.
3. Lubricate the distal tip of the Chartis Precision XL Catheter with a sterile lubricant for ease in passing the Chartis Precision XL Catheter through the working channel of the bronchoscope.

4. Depress plunger of sterile 3 mL syringe and reattach to one-way stopcock. Ensure stopcock is in open position.
5. Apply a vacuum to the balloon lumen by retracting the plunger of the 3 mL syringe and closing the one-way stopcock prior to inserting the catheter into the bronchoscope.
6. Carefully insert the Chartis Precision XL Catheter with the obturator into the working channel of the bronchoscope.
7. Position the bronchoscope into the target airway (e.g., RUL B1).
8. Advance the Chartis Precision XL Catheter approximately 2 cm beyond the distal tip of the bronchoscope until the catheter's distal balloon tip is visualized.
NOTE: Use caution to avoid kinking the catheter at the point where it is inserted into the bronchoscope valve.
9. Position the distal tip of the catheter into the target airway. Small positioning adjustments with the bronchoscope and/or the Chartis Precision XL Catheter should be made to ensure the balloon is positioned in the target airway location. The black proximal balloon marker identifies the proximal edge of the balloon.
10. Remove the obturator from the catheter.
11. Connect the Connector Set to the main lumen port on the catheter handle.
12. Flush the catheter.
 - Rotate the three-way stopcock such that the flush port is in the "OFF" position.
 - Retract the plunger of the 10 mL syringe and connect it to the flush port of the three-way stopcock. See Figure 3 below.
14. Confirm via visual inspection that the target airway is occluded. Air bubbles at the balloon should not be visible during respiration.
15. The system is now ready for data capture. See the User Manual for the Chartis Console for complete instructions for use.
16. Hold the bronchoscope/Chartis Precision XL Catheter in place while the Chartis Console takes measurements.
NOTE: If the Chartis Precision XL Catheter becomes blocked or clogged during use, insert the obturator to clear the lumen or small amounts of air may be injected through the inner lumen of the Catheter as described in the "Flush the catheter" section above.
17. When a measurement is complete (normally less than 5 minutes at each isolated lung segment), open the one-way stopcock and retract the plunger on the inflation syringe to deflate the balloon on the Chartis Precision XL Catheter. Use the syringe to pull additional vacuum as needed by detaching the syringe from the stopcock, purging air within the syringe, then reattaching to retract the plunger.
18. Re-insert the obturator.
19. Retract the Chartis Precision XL Catheter into the distal tip of the bronchoscope.
20. Reposition the bronchoscope near the next segment or lobe and repeat steps above to assess other airways as required.
21. When all measurements are complete, deflate the balloon, insert the obturator and remove the Chartis Precision XL Catheter by retracting it through the working channel of the bronchoscope.
CAUTION: To prevent balloon rupture, do not inflate the balloon more than 5 times for a single device.
22. Dispose of the Chartis Precision XL Catheter per institutional standard for biohazardous materials.
23. Any serious incidents or serious adverse events that occur due to the use of the Chartis Precision XL Catheter should be reported to Pulmonx (the manufacturer) and any relevant regulatory bodies (e.g., the European Competent Authority in the relevant EU Member State).


Figure 3: Catheter Flush Port Connection

- Depress the plunger of the 10 mL syringe to build up pressure.
 - Rotate the three-way stopcock valve to allow flow from the 10 mL syringe to the catheter. Do not direct flow to the console.
 - Detach the 10 mL syringe from the three-way stopcock and repeat the previous steps above to flush the catheter a second time.
 - Once flushing is complete, rotate the three-way stopcock valve such that the flush port is in the "OFF" position, ensuring free flow from the catheter to the console.
NOTE: Do not apply vacuum to the catheter or console during the flushing process.
13. Use the 3 mL syringe to inject up to 4 mL (cc) of air as needed for balloon expansion to create a seal in the airway being visualized and close the one-way stopcock to secure the inflated balloon. Do not over-inflate balloon in smaller airways.
CAUTION: To prevent balloon rupture, do not inject more than 4 mL (cc) of air.
- ## 8.0 HOW SUPPLIED
- The Chartis Precision XL Catheter is supplied sterile in a peel-open package.
 - The device is intended for single-patient use only. Do not re-sterilize the device.
 - Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken.
- ## 9.0 STORAGE AND HANDLING
- Storage: Store at room temperature in a clean and dry environment. Keep dry.
 - Protect Catheter from UV light exposure.
 - Handle with care.
 - Avoid extremes of temperature and humidity.
- ## 10.0 STERILE
- Sterilized with ethylene oxide. Do not use if the package is damaged or open.

Graphic Symbols Contained in Device Labeling

	Batch Code		Do Not Resterilize
	Catalog Number		Do Not Re-use
	Contents		Manufacturer
	Non-Pyrogenic		Date of Manufacture
	Sterilized Using Ethylene Oxide		Consult Instructions for Use
	Use-by Date		Keep Dry
	Medical Device		Unique Device Identifier
	Temperature Limit		Importer
	MR Unsafe		Do Not Use if Package is Damaged or Opened
	Authorized Representative in Switzerland/CH		Authorized Representative in the United Kingdom/UK
	Authorized Representative in the European Community/EU		Minimum Bronchoscope Working Channel Inner Diameter
	Caution, consult accompanying documents		Single Sterile Barrier System with Protective Packaging Outside
	Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.		
	Maximum Inflation Volume		

This product and/or its use are covered by one or more of the following United States patents: 7,883,471; 8,496,006; 8,523,782; 8,454,527; 8,808,194; 9,050,094; 9,107,606; 9,364,168; 9,533,116; 10,076,274; 10,314,516; 10,314,992; 10,413,244; 10,631,758; 10,758,239. Other U.S. patents pending.

This product and/or its use are covered by one or more of the following international patents: AT1435833; BE1435833; CH1435833; DE1435833; FR1435833; GB1435833; IE1435833; IT1435833; NL1435833; TR1435833; DE1838217; EP1838217B1; FR1838217; GB1838217; IE1838217; DE2285442; EP2285442B1; FR2285442; GB2285442; NL2285442; DE2614853; EP2614853B1; FR2614853; GB2614853; JP5430855B2; JP4301945; JP5452238B2. Other international patents pending.

Pulmonx Corporation
700 Chesapeake Drive
Redwood City, CA 94063
United States of America
Tel: +1 (650) 364-0400

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MDSS-UK RP Ltd.
6 Wilmslow Road, Rusholme
Manchester M14 5TP
United Kingdom

Pulmonx International Sàrl
Avenue de la Gare 2, 2000
Neuchâtel, Switzerland
Tel: +41 32 475 2070

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© 2025 Pulmonx Corporation 또는 그 계열사. All rights reserved. 모든 상표는 해당 소유자의 자산입니다.
Pulmonx GmbH