CHARTIS[®] PULMONARY ASSESSMENT SYSTEM CONSOLE

USER MANUAL

Manufactured by:

 Pulmonx Corporation
700 Chesapeake Drive
Redwood City, CA 94063
United States
Tel: +1 (650) 364-0400

European Authorized Representative:

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Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

CHARTIS[®] CONSOLE USER MANUAL

Copyright © 2023 by Pulmonx Corporation. All Rights Reserved.

This document may not be altered in any way without prior written consent of Pulmonx Corporation. Unauthorized reproduction, distribution, or use of this document, or any portion thereof, is strictly prohibited and may result in severe civil and criminal penalties.

Chartis[®] is a trademark of Pulmonx Corporation.

This product and/or its use are covered by one or more of the following United States patents: 7,883,471; 8,342,182; 8,454,527; 8,523,782; 8,808,194; 9,107,606; 9,364,168; and 9,439,583. Other U.S. patents pending.

This product and/or its use are covered by one or more of the following international patents: DE1435833B1; FR1435833B1; GB1435833B1; AT1435833B1; BE1435833B1; CH1435833B1; IE1435833B1; IT1435833B1; NL1435833B1; TR1435833B1; JP4301945B2; DE1838217B1; FR1838217B1; GB1838217B1; CH1838217B1; IE1838217B1; JP5430855B2; DE2056715B1; FR2056715B1; GB2056715B1. Other international patents pending.

Contents

1.	Introduction.....	5
2.	Indications for Use.....	5
3.	Contraindications.....	5
4.	Warnings and Precautions	6
5.	Hardware	7
6.	Security.....	8
7.	Symbols.....	9
7.1	Patient Interface Box (PIB)	11
7.2	Touchscreen Tablet Computer	12
8.	Getting Started	14
8.1	Console Preparation	14
8.2	Login Screen.....	16
8.3	Home Screen	17
8.3.1	Battery Charge Indicators	18
8.3.2	Language Selection.....	18
8.4	Setup Screen	19
8.4.1	Mode Selection	19
8.4.2	Date and Time Format.....	19
8.4.3	Warning Display.....	19
8.4.4	Real-time Airway Resistance display.....	20
8.4.5	Segmental Assessments	20
8.4.6	Default and Factory Settings	20
8.4.7	Setting Date and Time	20
8.4.8	System Information.....	21
9.	Patient Assessment.....	22
9.1	Patient Information.....	22
9.2	Select Airway Screen.....	23
9.2.1	Lobar Assessment.....	24
9.2.2	Segmental Assessment.....	24
9.3	Standard Mode Assessment Screen	26
9.4	Performing an assessment (Standard Mode)	29
9.4.1	Parameters	30
9.4.2	Window Buttons.....	30

9.4.3	Warnings.....	31
9.4.4	Rndx Trend.....	32
9.4.5	VT20 Trend.....	33
9.5	Post-Assessment Review (Standard Mode).....	34
9.6	Performing additional assessments.....	36
9.7	Ventilator Mode Assessment Screen.....	37
9.8	Performing an Assessment (Ventilator Mode).....	38
9.8.1	Parameters.....	38
9.8.2	VT20 Trend.....	39
9.9	Post-Assessment Review (Ventilator Mode).....	39
9.10	Patient Review Screen.....	40
9.10.1	Static Review.....	40
9.10.2	Dynamic Review.....	42
9.10.2.1	Video Recording.....	43
9.10.3	Video Review Screen.....	44
9.10.3.1	Video Player Screen.....	45
10.	Exporting Patient Assessment Data.....	46
11.	Interpreting Assessments.....	49
12.	Deleting Patient Assessment Data.....	57
13.	Console Shutdown.....	57
14.	Additional Information on the Tablet Computer.....	58
14.1	Power and Battery LED Status.....	58
14.2	Automatic Power On Feature.....	58
15.	Maintaining the Console.....	59
15.1	Preventive Maintenance.....	59
15.2	Cleaning Instructions.....	59
15.3	Storage.....	59
15.4	Caring for the Tablet Battery.....	59
15.5	Servicing.....	60
15.6	Internal Sensor Calibration.....	60
16.	Specifications.....	61
17.	Electromagnetic Compatibility (EMC) Compliance Information.....	62
18.	Troubleshooting.....	65
19.	Unit Conversion.....	68

1. Introduction

The Pulmonx Chartis Pulmonary Assessment System is designed to measure pressure and flow in order to detect the presence of collateral ventilation in isolated lung compartments. The System is designed to be used as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms. The System is to be operated only by trained bronchoscopists.

The system consists of:

- Chartis Catheter and Connector Set (see Chartis Catheter Instructions for Use)
- Chartis Console

The Chartis Console is a 2-part device consisting of a touchscreen tablet computer and a Patient Interface Box (PIB). The touchscreen tablet computer provides the user interface while the PIB houses the sensing equipment and associated electronics.

The Chartis Console includes:

- 1 Patient Interface Box (PIB)
- 1 Touchscreen tablet computer
- 2 Power Supplies with regional AC power cable
- 1 User Manual (this document)

Note: In this manual the word “device” is used to indicate the 2-part Chartis Console only, and the word “System” is used to indicate the Chartis Console in combination with the Chartis Catheter.

2. Indications for Use

The Chartis System is indicated for use by bronchoscopists during a diagnostic bronchoscopy in adult patients with Chronic Obstructive Pulmonary Disease (COPD) and emphysema in a bronchoscopy suite. The system, composed of the Chartis Catheter and Chartis Console, is designed to measure pressure and flow in order to calculate resistance to airflow and quantify collateral ventilation in isolated lung compartments. The Chartis Catheter is used through the working channel of a bronchoscope and connects to the Chartis Console. The Chartis Console is capital equipment that is reusable and displays the patient information.

3. Contraindications

The Chartis System is contraindicated in the presence of active infection or major bleeding diathesis. There are no known interfering substances.

4. Warnings and Precautions

- Check equipment for physical damage prior to use.
- As with any medical equipment, carefully route electrical cabling (i.e., power supply cabling) and pneumatic tubing (i.e., Chartis Catheter) to reduce the possibility of patient or user entanglement.
- To ensure patient safety, do not place the device in any position that might cause it to fall on the patient, or any position where its fall might cause it to pull on the Chartis Catheter while deployed in a patient.
- Use caution not to dislodge power cord when moving the Patient Interface Box while in use.
- Do not apply excessive pressure to the display screen.
- Do not subject the device to excessive shock or vibration.
- Do not subject the device to excessive magnetic fields. Damage to the touchscreen tablet computer and incorrect device operation may occur.
- Only the power supplies listed in the Specifications section of this manual are approved for use with the device. No other power supply should be used under any circumstances.
- The Chartis Catheter is the only type of catheter that can be used to safely interface with the Chartis Console.
- The Chartis Console needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this manual.
- Portable and mobile RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Chartis Console, including cables specified by Pulmonx. Otherwise, degradation of the performance of this equipment could result.
- Radio frequency transmitters such as Radio Frequency Identification (RFID) readers or systems at certain frequencies (e.g., 134.2 kHz and 13.56 MHz) potentially can cause electromagnetic interference (EMI) to the Chartis Console resulting in degradation of performance caused by an electromagnetic disturbance. Contact Pulmonx if any device malfunctions occur.
- The use of accessories, transducers and cables other than those specified by Pulmonx may result in increased electromagnetic emissions or decreased electromagnetic immunity of the Chartis Console and result in improper operation.
- No modification of this equipment is allowed.
- The Chartis Tablet Console is MR Unsafe. Keep away from magnetic resonance imaging (MRI) equipment.

5. Hardware

Warranty

Refer to the Chartis® Pulmonary Assessment System Limited Warranty card included with your shipment.

Disposal

For proper disposal of the Chartis Console or any of its accessories, dispose it in accordance with local, state, and federal regulations. Alternatively, you can return it to Pulmonx.

Radio Frequency Wireless Technology

The radio frequency functionality within the Chartis Console is completely redundant providing an alternative connection between the Patient Interface Box (PIB) and the touchscreen tablet computer which can be physically connected to each other via the docking station provided on the Patient Interface Box (PIB). The data transmitted wirelessly are not life-supporting or life-sustaining and do not contain any patient data. This device complies with part 15 of the FCC Rules. Operation is subject to the following conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: “Harmful interference” is defined by the FCC as follows,

Any emission, radiation or induction that endangers the functioning of a radio navigation service or of other safety services or seriously degrades, obstructs or repeatedly interrupts a radiocommunications service operating in accordance with FCC rules.

Below are the technical features of the wireless technology used:

Wireless Specification		
	Bluetooth Module in PIB	Bluetooth Module in Tablet
Design Name	BTM431 Laird Technologies	Intel AX210NGW
Transmitter Module FCC ID	PI4411B (Full Modular Approval)	YE3600-AX210NG
Bluetooth® Version	2.0+EDR	5.2
Transmit Class	Class 2	Class 2
Bluetooth SIG website Declaration ID	B016510	D043115
Frequency	2.402 – 2.480 GHz	2.402 – 2.480 GHz
Channels	79 channels Frequency Hopping, Adaptive Frequency Hopping	73, 68, 61 channels Frequency Hopping
Max Transmit Power	+0.0025 W (+4 dBm) from integrated antenna	+0.0093 W (+9.7 dBm)
Bluetooth Stack	V2.0 compliant. Fully integrated	V5.2 compliant
EU Directive	RED 2014/53/EU	RED 2014/53/EU

6. Security

Information Security

Patient information is private, confidential, and subject to various regulations. Access to the Chartis Console is designed to be controlled and limited to authorized users, and any patient data contained or generated by this device should be protected accordingly.

When entering data using the Chartis Console, it is your responsibility to protect your security credentials (e.g., passwords) and keep it confidential. Please do not use any patient protected health information (PHI) or individually identifiable health information of patients (e.g., names) to identify the assessments.

Wireless Security

The Chartis Patient Interface Box (PIB) and the Chartis touchscreen tablet computer are directly connected when the tablet computer is placed in the docking station of the Patient Interface Box. The Chartis touchscreen tablet computer may also be operated 'undocked', allowing it to connect with the Patient Interface Box via Bluetooth[®] Version 2.0+EDR using frequencies within 2.402 – 2.480 GHz - Class 2. To hinder connections with untrusted Bluetooth device(s) in the vicinity and to prevent malicious parties from intercepting the Bluetooth signal, multiple layers of security have been incorporated:

- The Chartis Tablet is not a discoverable Bluetooth device.
- The Chartis Tablet searches for the PIB using a unique prefix.
- Once discovered, the Chartis Tablet pairs with the PIB using a preset PIN number.

Confidentiality

Confidentiality of patient information is assured as follows:

- The data transmitted between the Chartis Tablet and the Chartis PIB does not contain any patient-identifying information.
- All patient assessments are stored within the Chartis Tablet and access is limited to specific user login credentials.
- In the case of multiple users, you cannot access assessments created under another user's login credentials.
- Exported reports are also password protected and require the same password that was used to login to the user account when the patient report was exported.

Wireless Data Integrity

The Chartis software incorporates data integrity checks for data transmitted/received via the Bluetooth to ensure the completion and validity of the data received. If any data is incomplete or invalid, an error is displayed.

7. Symbols

In this manual, displayed on device labels or the device itself, the following symbols are used:

Symbol	Meaning
	Contents
	Catalog Number
	Model Number
	Serial Number
	Date of Manufacture
	Manufacturer
	European Authorized Representative
	UK Authorized Representative
	Swiss Authorized Representative
	Importer
	BF (Body Protected). The input connectors are suitable for connection to humans provided that there is no direct connection to the heart.
	See accompanying documentation. The supplied documentation must be consulted for operating, cautionary or safety information before using the device.
	Caution

Symbol	Meaning
	Non-sterile
	WEEE (Waste from Electrical and Electronic Equipment). Device must be taken to a recycling center for proper disposal.
	Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.
	Fragile, handle with care
	Keep Dry
	Temperature Limit
	Humidity Limitation
	This way up
	Stacking Limitation: Do not stack more than two high when placing on a pallet.
	To indicate that the marked item or its material is part of a recovery or recycling process.
	MR Unsafe
	Medical Device

7.1 Patient Interface Box (PIB)

Symbol	Meaning
	“ON”/”OFF” (push-push). Push to turn power ON or OFF
	Model Number
	Serial Number
	BF (Body Protected). The input connectors are suitable for connection to humans provided that there is no direct connection to the heart.
	Direct Current. The device uses direct currents.
	Date of Manufacture
	WEEE (Waste from Electrical and Electronic Equipment). Device must be taken to a recycling center for proper disposal.
	Non-sterile
	Manufacturer
	European Authorized Representative
	UK Authorized Representative
	Swiss Authorized Representative
	Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.
	See accompanying documentation. The supplied documentation must be consulted for operating, cautionary or safety information before using the device.

Symbol	Meaning
	MR Unsafe
	Medical Device
	Use By Date

7.2 Touchscreen Tablet Computer

Symbol	Meaning
	“ON”/“OFF” (push-push). Push to turn power ON or OFF
	Model Number
	Serial Number
	Direct Current. The device uses direct currents.
	Date of Manufacture
	WEEE (Waste from Electrical and Electronic Equipment). Device must be taken to a recycling center for proper disposal.
	Manufacturer
	European Authorized Representative
	UK Authorized Representative
	Swiss Authorized Representative

Symbol	Meaning
	<p>Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.</p>
	<p>Non-sterile</p>
	<p>USB Type-A</p>
	<p>USB Type-C</p>
	<p>See accompanying documentation. The supplied documentation must be consulted for operating, cautionary or safety information before using the device.</p>
	<p>MR Unsafe</p>
	<p>Medical Device</p>

8. Getting Started

8.1 Console Preparation

Read all instructions carefully and in their entirety.

The Chartis Console may be placed on a tray, desktop or cart.

Warning No device is approved for connection to the touchscreen tablet computer 3.5mm audio I/O jack, LAN port, and smart card reader.

TABLET COMPUTER SETUP FOR THE FIRST TIME:

To power on the touchscreen tablet computer for the first time, connect it to power either via docking it into the powered ON Patient Interface Box or by connecting it directly to the power outlet using the provided power supply.

Note: The tablet computer is set into a hibernation mode during shipping to prevent the tablet battery from draining. Thus, out of the box, the tablet computer is not responsive to the pushing of the power button on the side of the tablet. The tablet's hibernation mode is disabled after the first connection to power either via docking into the powered on PIB or connecting the tablet to the power outlet via the power supply. After the tablet's hibernation mode is disabled, the tablet computer's power button may be used as one method to power on the tablet.

CHARTIS CONSOLE SET UP CONFIGURATIONS

SETUP CONFIGURATION 1 – Tablet docked on PIB

1. Connect the Patient Interface Box (PIB) to a power outlet using the power supply provided.
2. Connect the touchscreen tablet computer to the PIB by docking the tablet computer into the docking station on the PIB.
3. Turn ON the system by pressing the power button on the front panel of the PIB. Both the PIB and touchscreen tablet will be powered on.
4. The tablet software will display a data connection symbol  to indicate that a data connection was established between the PIB and the touchscreen tablet computer.

SETUP CONFIGURATION 2 – Undocked Bluetooth

Note: Setup Configuration 2 may only be used when the tablet battery is sufficiently charged.

1. Connect the Patient Interface Box (PIB) to a power outlet using the power supply provided.
2. If connected, disconnect the touchscreen tablet computer from the PIB by removing the tablet computer from the docking station on the PIB.
3. Turn ON the touchscreen tablet computer using the power button  on the side of the tablet computer.
4. Turn ON the PIB by pressing the power button on the front panel.
5. The tablet software will display a data connection symbol  to indicate that a data connection via Bluetooth was established between the PIB and the touchscreen tablet computer.

SETUP CONFIGURATION 3 – Tablet docked using battery power

Note: Setup Configuration 3 may only be used when the tablet battery is sufficiently charged.

1. Connect the touchscreen tablet computer to the Patient Interface Box (PIB) by docking the tablet computer into the docking station on the PIB.
2. Turn ON the Touchscreen tablet computer using the power button  on the side of the tablet computer.
3. Turn ON the PIB by pressing the power button on the front panel.
4. The tablet software will display a data connection symbol  to indicate that a data connection was established between the PIB and the touchscreen tablet computer.

The Chartis Catheter will connect to the front quick-disconnect socket via the Connector Set. The socket on the back of the console is to remain unplugged during console operation.

Warning *The Chartis Catheter should only be connected to the front quick-disconnect socket of the Patient Interface Box.**


Precaution *All dust plugs must be unplugged before initiating a patient assessment.*

* Refer to the Chartis Catheter Instructions For Use for preparation and proper attachment of the catheter to the console via the connector set.

CHARTIS CONSOLE DATE AND TIME SET UP


During the first power-up of the Chartis Console, the following screen will be displayed. Using the on-screen month, day, year, hour, minute, AM/PM buttons, set the system clock to the current date and time and press the **OK** button.

Additionally, if the Chartis Console displays this screen at power-up, you must set the system clock to the current date and time before performing a patient assessment to ensure that patient records reflect the correct date and time data.


8.2 Login Screen


Once the system software finishes booting up, the LOGIN screen will appear. On the LOGIN screen, an existing user may login by entering an existing username and password and then pressing the *Login* button.

Additionally, a new user may be created by pressing the *New User* button and creating a new username and password.


If the Chartis Console displays the following message after logging into a user account, you must reset the date and time before performing a patient assessment to ensure that patient records reflect the correct date and time data.

Dismiss this message by pressing **OK**, and then follow the instructions in Section 8.4.7 of the User Manual for setting the correct date and time.

8.3 Home Screen

Once a User is logged in, the HOME screen will appear.

8.3.1 Battery Charge Indicators

The status of the battery charge is represented by the battery symbol:

- When the charge is in the 40-100% range, the color is white: 
- When the charge is in the 20-39% range, the color is orange: 
- When the charge is in the 0-19% range, the color is red: 

Next to the battery symbol is a percentage value of the current charge. When the tablet battery is being charged from an external power supply (PIB dock or standalone tablet charger), the lightning bolt symbol is displayed: 


See section 14 (Additional Information on the Tablet Computer) for definition of LED status icons located on the left edge of the tablet computer screen.

8.3.2 Language Selection


The default language for the Chartis Console is English. To change the displayed language on the Console, follow the directions below.


From any screen, touch the picture of the country flag. A language selection screen will pop up. Select the country flag for the desired language. This change will take place instantly and will remain in effect even when the console is power cycled.

Note: The additional language screens for data entry will use the English keyboard.

8.4 Setup Screen

From the HOME screen, press the  icon on the bottom left corner of the page to display the SETUP screen. The Service Mode screen (not shown) is only accessible via correct password entry and available to Chartis Console manufacturing representatives for calibration, testing, and troubleshooting.

8.4.1 Mode Selection

The MODE SELECTION checkboxes allow you to choose between Standard Mode and Ventilator Mode [default] for assessing a patient.

8.4.2 Date and Time Format

The DATE check box allows you to choose the date format that is displayed throughout the interfaces on the console. You can choose between DD/MM/YYYY (default) and MM/DD/YYYY. The TIME check box allows you to choose the time format that is displayed throughout the interfaces on the console. The selection options are in the 12 hour or 24 hour formats.

8.4.3 Warning Display

The “Display Warnings” checkbox allows you to enable/disable warnings during a Standard Mode assessment, as well as on Review screens and in assessment reports. You can turn on the warnings by checking the box. To view a list of warnings, see Section 9.4.3.

8.4.4 Real-time Airway Resistance display

The “Display Rrt graph” checkbox allows you to enable/disable the real-time airway resistance display during a Standard Mode assessment, as well as on Review screens and in assessment reports.

8.4.5 Segmental Assessments


The “Show Segmental” checkbox allows you to enable/disable the segmental lung map on the assessment lung map.

8.4.6 Default and Factory Settings

The *Save as Default* button allows you to save your settings so that they are the same when you restart the console. Choose your settings, and then press the *Save as Default* button to save them. If you would like to revert back to the factory settings, press the *Factory Reset* button.


8.4.7 Setting Date and Time

Press *Set Date/Time* to open the Set Date and Time screen. Press (+) or (-) buttons on each value to set the calendar date and local time. When the correct date and time are set, press the *OK* button.

8.4.8 System Information


From the **SETUP** screen, press the *System Information* button to display the system information (sample screen shown below). Press *Close* to return to the **HOME** screen.

9. Patient Assessment


9.1 Patient Information

From the HOME screen, press **Assess** to display the NEW PATIENT screen below. Press the field next to Patient ID to enter the Patient ID using the keyboard below.


Once the Patient ID has been entered, press the **Assess** button under the Patient ID textbox on the screen to create the patient record and continue to the **SELECT AIRWAY** screen.

Note: Patients previously assessed that already have patient records can be reassessed through the **PATIENT REVIEW** screen. From the HOME screen, press **Review**. Press a row to highlight the patient you would like to do another assessment for, and press **Assess**.

9.2 Select Airway Screen

The **SELECT AIRWAY** screen for lobar and segmental assessments appears after pressing **Assess** on the **NEW PATIENT** screen. *Note: the Segmental Airway Selection will only be available if this setting is checked on the SETUP screen.*


The **MODE SELECTION** check boxes toggle between *Ventilator* and *Standard* modes. The **AIRWAY** selection check boxes toggle between the *Lobar* and *Segmental* lung maps. Selecting *Lobar* will record assessments at the lobe level and *Segmental* will allow assessments of individual airways within the lobe. The number of assessments in each mode performed for each lobe and segmental region is shown below the lobe name with the following format:

[# LOBAR ASSESSMENTS]

[# SEGMENTAL ASSESSMENTS]


9.2.1 Lobar Assessment

Press *Lobar* using the AIRWAY selection check boxes [default selection]. Select the mode to perform the assessment in (Standard or Ventilator). Select a lobe to assess by pressing the desired lobe in the diagram. The ASSESS screen will appear.


9.2.2 Segmental Assessment

Press *Segmental* using the AIRWAY selection check boxes (as shown below).

Select a lobe to view the segmental lung map for that lobe by pressing on the desired lobe in the diagram. The system will navigate to the segmental lung map (as shown below). Select the mode to perform the assessment in (Standard or Ventilator). Select a segment to assess by pressing on the desired segment in the diagram. The ASSESS screen will appear.


9.3 Standard Mode Assessment Screen

The selected lobe and corresponding total number of assessments performed for that lobe will be displayed at the top left of the assessment screen in the format <LOBE> <ASSESSMENT #>. The Standard Mode ASSESS screen has three charts: the upper chart displays flow and pressure data, while the lower charts display the Resistance data or VT20 data depending on which is selected. To view the two different lower charts, the user may toggle between the two graphs, marked by the red button located between the upper and lower charts.


CHART	PARAMETER / COLOR	DESCRIPTION
UPPER	F (mL/min) / <i>ORANGE</i>	EXPIRATORY FLOW (F) VS. TIME
	P (cmH ₂ O) / <i>LIGHT BLUE</i>	INSPIRATORY PRESSURE (P) VS. TIME
LOWER	Rndx / <i>GREEN</i>	AIRWAY RESISTANCE INDEX (Rndx) VS. TIME (UNITS OF cmH ₂ O×sec/mL)
	Rrt / <i>GREEN (must be enabled in settings screen)</i>	REAL-TIME AIRWAY RESISTANCE (UNITS OF cmH ₂ O×sec/mL)
LOWER	VT20 (mL) / <i>YELLOW</i>	VOLUME TREND FOR PREVIOUS 20 SECONDS (VT20)


Display of the Resistance chart:

Display of the VT20 Chart:

The baseline state of the combined flow and pressure trace is displayed below. When the dataset reaches the far right of the screen, it continues from the left, overwriting the previous data sweep. Note that the previous data sweep is dimmed to highlight the location of the newest data (grey vertical cursor).

9.4 Performing an assessment (Standard Mode)

Press the *Start* button on the lower left hand corner to initialize the assessment. This collects 10 seconds of baseline flow and pressure data from the isolated airway.

After the 10-second initialization is complete, the console will close a valve to prevent flow back into the isolated lobe during inspiration. This will result in the inspiratory pressure increasing in magnitude.

- The top chart continuously auto scales the positive and negative regions of the vertical axis based on the maximum value of flow and inspiratory pressure, respectively, displayed throughout the assessment.
- The resistance chart on the bottom has a vertical scale fixed at a maximum Rndx of 20 cmH₂O×sec/mL.
- The VT20 chart on the bottom continuously auto scales the vertical axis based on the maximum value of VT20 displayed throughout the assessment.

9.4.1 Parameters

At the left side of the ASSESS screen, the following parameters are displayed in Standard Mode:


PARAMETER	DESCRIPTION	UNITS
Time	ASSESSMENT DURATION	minutes:seconds
Flow (mL/min)	POSITIVE FLOW	mL/min
Flow Trend	EXPIRATORY FLOW VOLUME AS A PERCENTAGE OF THE AVERAGE EXPIRED VOLUME PER BREATH CALCULATED DURING THE 10-SECOND INITIALIZATION PERIOD	%
VOL (mL)	TOTAL POSITIVE AIR FLOW VOLUME SINCE ASSESSMENT START	mL
VT20	CURRENT VOLUME EXPIRED MINUS THE TOTAL VOLUME EXPIRED FROM TWENTY SECONDS PREVIOUS	mL
RNDX	AIRWAY RESISTANCE INDEX (Rndx) VS. TIME	cmH ₂ O×sec/mL

9.4.2 Window Buttons

BUTTON	DESCRIPTION
	The ZOOM button, located between the upper and lower charts, allows the user to zoom in (+) or zoom out (-).
	Press the START button to start an assessment.
	Press the STOP button to stop at the end of a satisfactory assessment period. The REVIEW screen (post-assessment) will appear.
	Press the Reset Rndx button to reset the resistance index calculation. A caret (^) shall be displayed in the pressure and Rndx graphs, indicating when the user reset the calculation.

9.4.3 Warnings

During an assessment, you may see the following warnings flash on the display next to the icon between the top and bottom charts (see below) if “*Display Warnings in Assessments and Reports (Standard Mode)*” is enabled in the SETUP screen.

The following table describes warnings that may be encountered during a Standard Mode assessment.

WARNING	DESCRIPTION
Clog Detected	The system displays a text warning and red warning indicator bars at the bottom of the chart to indicate detection of clogs during a Standard Mode assessment.
No Tidal Breathing Detected	The system displays a text warning and red warning indicator bars at the bottom of the chart to indicate failure to detect tidal breathing during the 10-second pre-assessment and post-assessment periods of a Standard Mode assessment.

If enabled on the SETUP screen, the warning bars will be displayed on the Review screens and assessment reports.

9.4.4 Rndx Trend

The 2-image sequence below shows the rise in calculated airway resistance index (R_{ndx}) that may occur during a Standard Mode assessment.


The **Stop** button is pressed to end an assessment, and the valve then opens during inspiration, as illustrated by the drop in magnitude of the inspiratory pressure. Ten seconds of post-assessment flow and pressure data are recorded.

9.4.5 VT20 Trend

The image below shows the declining VT20 line that may occur during a CV-, Standard Mode assessment.

9.5 Post-Assessment Review (Standard Mode)


The REVIEW screen (post-assessment) appears showing the pre-assessment initialization, assessment period, and post-assessment period, all separated by two vertical grey dashed lines.


Pressure and Flow data (upper chart): Touch and drag the grey cursor across the screen to determine the value of flow (orange) or pressure (light blue) anywhere in the patient assessment dataset.

Resistance Index (lower chart): Shows the time history of the airway resistance index (green) and real-time resistance (light blue). Touch and drag the grey cursor across the screen to determine the value of airway resistance index, real-time airway resistance (if enabled), and VT20 anywhere in the patient assessment dataset.

Assess Time: Reports the duration of the assessment.

Vol (mL): Reports the total volume exhaled during the assessment.


If there are reasons to discard the assessment, press *Discard* and then press *Yes* to confirm discard.

Otherwise, press *Accept* to accept the assessment.

The system will redirect to the most recently visited airway selection screen (lobar or segmental map). The user can select a new airway for assessment, or the user can press the *Home* button to return to the HOME screen (return to Section 8.3).

9.6 Performing additional assessments

If necessary, deflate the balloon, reposition the catheter, and re-inflate per the Chartis Catheter Instructions for Use. Upon returning to the SELECT AIRWAY screen, the number of assessments will update accordingly.

Perform additional assessment(s) per Section 9.3, STANDARD MODE ASSESSMENT SCREEN.

9.7 Ventilator Mode Assessment Screen

Ventilator Mode is a simplified assessment screen to display flow and VT20 data when assessing patients under a ventilator. Presence of continuous flow determines presence of CV. As a result, inspiratory pressure is no longer required on display. The valve actuation algorithm does not change between Ventilator and Standard modes.

Note: The Ventilator Mode provides a display option that may be preferred when the physician selects to perform the bronchoscopy and Chartis assessment on a patient under general anesthesia with mechanical ventilation. The term “Ventilator Mode” does not signify it is intended for use in patients that are ventilator dependent outside of the procedure.

The Ventilator Mode ASSESSMENT screen has two charts: the upper chart displays flow data and the lower chart displays VT20.

CHART	PARAMETER / COLOR	DESCRIPTION
UPPER	Flow (mL/min) / <i>ORANGE</i>	EXPIRATORY FLOW (F) VS. TIME
LOWER	VT20 (mL) / <i>YELLOW</i>	VOLUME TREND FOR PREVIOUS 20 SECONDS (VT20)

9.8 Performing an Assessment (Ventilator Mode)

The steps leading up to an assessment such as creating a patient (refer to Section 9.1 for Patient Information) and selecting a lobe (see Section 9.2 for Select Airway Screen) are the same as Standard Mode. The assessment screen for Ventilator Mode is similar to the Standard Mode except it does not show the combined Flow and Pressure trace (see 9.3). The upper chart shows only positive expiratory flow while the lower chart shows the VT20 (Volume trend for previous 20 seconds). The Peak Flow Trend that is displayed in Ventilator Mode is the Maximum Flow in every second of displayed flow signals expressed as a percentage of the average peak flow calculated in the 10-second pre-assessment period. Navigating to the REVIEW screen (post-assessment) is the same as for Standard Mode assessments (refer to Section 9.5).


9.8.1 Parameters

At the left side of the ASSESS screen, the following parameters are displayed in Ventilator Mode:

PARAMETER	DESCRIPTION	UNITS
Time	ASSESSMENT DURATION	minutes:seconds
Flow (mL/min)	POSITIVE FLOW	mL/min
Peak Flow Trend	MAXIMUM EXPIRATORY FLOW IN EVERY SECOND OF DISPLAYED FLOW SIGNALS	%
VOL (mL)	TOTAL POISTIVE AIR FLOW VOLUME SINCE ASSESSMENT START	mL
VT20	CURRENT VOLUME EXPIRED MINUS THE TOTAL VOLUME EXPIRED FROM TWENTY SECONDS PREVIOUS	mL


9.8.2 VT20 Trend

The image below shows the declining VT20 line that may occur during a CV-, Ventilator Mode assessment:

9.9 Post-Assessment Review (Ventilator Mode)

Follow the same steps for Standard Mode to review a Ventilator Mode assessment (refer to Section 9.5). The REVIEW screen (post-assessment) for Ventilator Mode will only display the summary Flow chart, VT20 chart, and Volume measurement.


9.10 Patient Review Screen

From the HOME screen, press *Review* to display the PATIENT REVIEW screen. The user may click *Assess*, *Static Review*, *Dynamic Review*, or *Video* (when a recording is made).


9.10.1 Static Review

Select a row to highlight the patient for review, and press *Static review*.


The **SELECT AIRWAY** screen appears. Select an airway to review by pressing on the desired area in the diagram. The **AIRWAY** check boxes can be used to review segmental assessments by pressing *Segmental*, pressing a lobe, and then pressing a segment to review (refer to Section 9.2.2).

The **STATIC REVIEW** screen appears, and the arrow buttons at the bottom of the screen can be used to toggle through multiple assessments for the same lobe. Touch and drag the grey cursors to determine the value of flow, pressure, airway resistance index, real-time airway resistance (if enabled), and VT20 anywhere in the patient assessment dataset.


In the sample image above, the assessment was **ACCEPTED**. In the sample image below, the assessment was **DISCARDED**.


To zoom in at the cursor location, press the magnifying glass (+) located between the two charts. Press the magnifying glass (-) to zoom out.

9.10.2 Dynamic Review

From the Patient Review screen, select *Dynamic review*. The user will then be directed to the SELECT AIRWAY screen. The user can select the airway and view a live playback of the assessment. The DYNAMIC REVIEW screen appears, and the arrow buttons at the bottom of the screen can be used to toggle through multiple assessments for the same lobe.


After the live assessment has been completed, the notification below will appear, indicating that the dynamic review of the assessment has finished.

9.10.2.1 Video Recording


In the DYNAMIC REVIEW screen, the user can select the  button, which will begin recording the assessment.

Once an assessment recording starts, a red box will appear on the screen, indicating that a recording is in progress. You may record from the start of an assessment until any desired stopping point. The zoom controls during the recording may be used to display various screens during the recording. Once the user wishes to end a recording, click


to stop the recording.


9.10.3 Video Review Screen

From the PATIENT REVIEW Screen, the user may select the *Video* button to view any video recordings they have recorded previously.


After selecting *Video* from the PATIENT REVIEW screen, the VIDEO REVIEW screen will appear. Here the user may play or delete the video recording(s).

9.10.3.1 Video Player Screen

If the user clicks *Play*, the VIDEO PLAYER screen will appear and the user may play, pause, or use the slider to playback different parts of the assessment recording.

10. Exporting Patient Assessment Data

Insert a USB flash memory device into the USB Type-A port or USB Type-C port on the side of the touchscreen tablet computer.

Note: Connect only one storage device to the tablet computer at a time.

Warning *No other type of device is approved for connection to the touchscreen tablet computer USB ports including, but not limited to, keyboards, mice or other pointing devices, powered memory devices (which require external power to operate), still cameras, video cameras, music players, cellular telephones, and other multimedia devices.*

From the PATIENT REVIEW screen, press the checkbox next to the patient (or anywhere on the row) to select it, and then press *Export*.

The patient data files, including a report in PDF format for each saved assessment (see below) and any video recordings made will be copied to the USB flash memory device.

PDF Report Filename Format:

[Patient ID]_[Lobe]-[Assessment #] [Assessment Mode] Report.pdf

Video Filename Format:


[Patient ID]_[Lobe]-[Assessment #] [Assessment Mode] Video.mkv

(Note: to open the video file, which is exported as an MKV file format, you must have H.265/HEVC codec installed or you will need to convert the file to a readable format.)

PDF Report and Video file Drive Location:

[Drive]:\Pulmonx\ChartisConsole\patient_data\[Patient ID]_[Assessment Mode]
[Assessment Mode] will be either STD (Standard Mode) or VEN (Ventilator Mode).

When exporting patient data from the Chartis Console, you must complete the data export process by pressing *OK* to this message **before** removing the USB flash memory device.


Warning *Removing the USB device before the data export process is complete will corrupt the data on the USB flash memory device.*

The reports can be opened on any Windows PC system with a PDF viewer. The report and zipped video file are password protected and require the same password that was used to login to the user account when the patient was exported.

Sample PDF reports for Standard Mode and Ventilator Mode are shown below:

SAMPLE STANDARD MODE ASSESSMENT REPORT

CHARTIS[®] CONSOLE USER MANUAL

SAMPLE VENTILATOR MODE ASSESSMENT REPORT

Ventilator Assessment Report



S/N 506596; B-S/N 506596; SW VER 6.1.4

Patient Information

Patient ID: TEST-PATIENT
Date of Assessment: 04/12/2023

Assessment

Assessed Airway: RML (Assessment # 1)
Assessment Start Time: 3:54 PM
Assessment Duration: 01:54
Total Exhaled Volume: 23 mL
Assessment Result: **Accepted**

11. Interpreting Assessments

An example output for an acceptable Chartis System assessment is shown below. While reviewing collateral ventilation status of an assessment report, there are several important items to note, including:

- 1) The report should demonstrate consistent tidal breathing prior to start of assessment. This indicates that the patient was adequately sedated. Ensuring adequate sedation will aid in achieving complete sealing of the assessed airway with the Chartis system.
- 2) Duration of the assessment should be 2 minutes or longer.
- 3) Throughout the assessment, consistent deeper negative pressure should be observed. This indicates that proper catheter balloon positioning and sealing was achieved.
- 4) At the end of the assessment, the total volume of air exhaled by the patient during the assessment should have been greater than 50 ml.
- 5) During the post-assessment period, ensure that tidal breathing pattern returned.
- 6) If the report indicates that the patient coughed, the catheter tip clogged, or the balloon seal was lost in such a manner that affect the assessment, then the assessment should be repeated.
- 7) During the assessment, if there was an instantaneous loss or drop in flow within the first 30 seconds and flow did not resume, the assessment should be repeated.


(Note: The figure below is NOT displayed on the Chartis application, all three charts are provided below to view Resistance chart alongside VT20 to help understand how to use these charts to determine if an assessment is CV- or CV+)

Standard Assessment of an airway that has Little or No Collateral Ventilation


1. When high resistance to collateral ventilation or high collateral resistance is noted, the Chartis assessment indicates that the assessed airway has little or no collateral ventilation (is CV negative or CV-).
2. An example of a report for a CV- assessment is shown below.
3. The key features of a CV- assessment include:
 - a. The peaks of the flow indicator waveforms (orange) show a steady downward trend.
 - b. The Airway Resistance index (Rndx) and / or Real-time Airway Resistance (Rrt) profile lines climb in a stepwise and steady fashion.
 - c. The VT20 has a steady downward trend toward 0.
 - d. Potential artifacts to avoid include:
 - i. Instantaneous climbing of the Rndx and Rrt lines should not be considered a CV- assessment. This observation could be an artifact of an underlying issue such as clogging of the catheter tip.
 - ii. An instantaneous drop in flow seen within the first 30 seconds of the assessment should not be considered a CV- assessment. This observation could be an artifact of distal airway collapse.

(Note: The figure below is NOT displayed on the Chartis application, all three charts are provided below to view Resistance chart alongside VT20 to help understand how to use these charts to determine if an assessment is CV- or CV+)

Ventilator Assessment of an airway that has Little or No Collateral Ventilation

1. In Ventilator Mode, if there a downward trend in flow, collateral ventilation is noted and the Chartis assessment indicates that the assessed airway has low collateral ventilation (CV negative or CV-).
2. An example of a report for a CV- assessment is shown below.
3. The key features of a CV- assessment include:
 - a. The peaks of the flow indicator waveforms (orange) show a steady downward trend.
 - b. The VT20 has a steady downward trend toward 0.
 - c. Potential artifacts to avoid include:
 - i. Instantaneous decrease of the VT20 lines should not be considered a CV- assessment. Do not reference the VT20 plot if the line turns dashed, this indicates that the VT20 line is no longer valid. The VT20 line could change to dashed pattern if there are any underlying issue such as clogging of the catheter tip.
 - ii. An instantaneous drop in flow seen within the first 30 seconds of the assessment should not be considered a CV- assessment. This observation could be an artifact of distal airway collapse.

Standard Mode Assessment of an airway that demonstrates presence of Collateral Ventilation

1. When low resistance to collateral ventilation or low collateral resistance is noted, the Chartis assessment indicates that the assessed airway has high collateral ventilation (CV positive or CV+).
2. An example of a CV+ assessment is shown below.
3. The key features of a CV+ assessment include:
 - a. The peaks of the flow indicator waveforms (orange) remain flat and do not show a downward trend.
 - b. The Airway Resistance index (Rndx) and / or Real-time Airway Resistance (Rrt) profile lines have little or no change and do not steadily increase.
 - c. The VT20 has little or no change and does not trend toward 0.


(Note: The figure below is NOT displayed on the Chartis application, all three charts are provided below to view Resistance chart alongside VT20 to help understand how to use these charts to determine if an assessment is CV- or CV+)

Other measurement anomalies which might affect the outcome of the collateral ventilation assessment include clogged catheter tip, instantaneous drop in flow, loss of balloon seal, and patient coughing. It should be noted that these anomalies can occur together or independently during the course of an assessment.

Ventilator Mode Assessment of an airway that demonstrates presence of Collateral Ventilation


1. In Ventilator Mode, if there is no downward trend in flow, collateral ventilation is noted and the Chartis assessment indicates that the assessed airway has high collateral ventilation (CV positive or CV+).
2. An example of a CV+ assessment is shown below.
3. The key features of a CV+ assessment include:
 - a. The peaks of the flow indicator waveforms (orange) remain flat and do not show a downward trend.
 - b. The VT20 line has little or no change and does not trend toward 0.

Airway assessment that has mucous clogging.


An example of a report with clogging at the catheter tip is shown below.

(Note: The figure below is NOT displayed on the Chartis application, all three charts are provided below to view Resistance chart alongside VT20 to help understand how to use these charts to determine if an assessment is CV- or CV+)


Airway assessment with instantaneous drop in flow.

An example of an assessment with an instantaneous drop in flow is shown below.


Airway assessment with loss of Chartis Catheter balloon seal.

An example of a report in which the balloon seal was compromised is shown below.


Airway assessment with inadequate Chartis Catheter balloon seal and coughing.

An example of a report with inadequate balloon seal and a cough during the assessment is shown below.

12. Deleting Patient Assessment Data

From the **PATIENT REVIEW** screen, press the checkbox next to the patient (or anywhere in the row) to select it, and then press *Delete*. Confirm that is the desired action to complete the deletion process.

13. Console Shutdown

To power down the console from the **HOME** screen, press the  button located on the lower right corner. Then, press the *Shutdown* button if ready to shutdown. Disconnect the Chartis Catheter from the front quick-disconnect socket. Insert the 2 dust plugs into their respective quick-disconnect sockets.

Precaution *When the device is not in use, the two (2) dust plugs should be inserted into their respective quick-disconnect sockets.*

Wait until the screen turns black, then push the switch on the Patient Interface Box (PIB) to complete the shutdown process. Turn the (PIB) Power Switch off within 30 seconds, whether or not the black screen appears.

The power supply must be unplugged from the wall power outlet in order to fully isolate the system from electrical power.

14. Additional Information on the Tablet Computer

14.1 Power and Battery LED Status

There are two LED status icons on the left edge of the tablet computer screen.

 LED indicates the tablet battery status, and the  LED indicates the tablet computer power status (ON or OFF).

The table below defines the power and battery LED statuses:

Icon	Power condition	LED Status	Definition
	Battery powered	Solid Blue Light	Battery charge level is between 26% to 100%
		Solid Orange Light	Battery charge level is between 10% to 25%
		Blinking Orange Light	Battery charge level is below 10%
	Tablet computer is connected to external power source (docked on powered on PIB or power supply is plugged into tablet)	Solid Blue Light	Battery charge level is 100% (fully charged)
		Blinking Blue Light	Battery is charging
		Fast Blinking Orange Light	Tablet battery temperature alarm
		Solid Brown Light (Both Blue and Orange LED are ON)	Tablet battery is dead
	N/A	Solid Blue Light	Tablet is ON
		No Light	Tablet is OFF

14.2 Automatic Power On Feature

The touchscreen tablet computer will power ON automatically under two conditions:

- Tablet computer is docked on a powered on Patient Interface Box (power button on the front panel of PIB is pressed inward).
- The power supply provided with the system is plugged directly into the tablet.

However, if the tablet computer is power cycled, wait 10+ seconds after tablet computer is OFF before attempting to power ON the tablet computer via the automatic power on methods described above. The 10+ second wait does not apply if powering ON the tablet computer using the power button on the side of the tablet computer.

15. Maintaining the Console

This section describes maintenance, storage, and service for the Chartis Console.

15.1 Preventive Maintenance

The Chartis Console may be surface-cleaned after each patient use.

Precaution *Clean the Chartis Console according to the instructions in this document.*

Precaution *Insert dust plugs into the quick-disconnect sockets before cleaning.*

15.2 Cleaning Instructions

Verify that the console is powered off. Use a lint-free cloth to clean the touch screen tablet computer and the Patient Interface Box. Spray a 70% isopropyl alcohol (IPA) solution onto the lint-free cloth, and then wipe down both the tablet computer and Patient Interface Box. Never spray cleaning solutions directly onto the device surfaces.

Precaution *Paper towels or other cloths that are not lint-free may leave a residue and/or scratch the touchscreen.*

Precaution *Glass cleaners that contain ammonia or alcohol can degrade the touchscreen. Tap water or mineral water could leave white marks on the touchscreen due to the dissolved salts.*

Precaution *Do not oversoak the lint-free cloth and do not spray the 70% IPA directly onto the patient interface box or touchscreen tablet computer to help minimize any cleaner from leaking into the console.*

15.3 Storage

When not in use, store the Chartis Console in a cool, dry environment. See Section 16 for temperature and humidity storage specifications.

15.4 Caring for the Tablet Battery

- It is recommended to keep the tablet battery charge level above 20%.
- Charge the tablet battery between uses.
- To prevent the tablet battery from swelling, avoid high heat and avoid letting it stay in a low-capacity state (battery charge level less than 20%) for too long.
- To maximize the performance of the tablet battery, condition the battery once a month. To condition the battery, run the tablet on battery power until the battery's charge level is below 20%. Then, use AC power until the battery is fully charged. Avoid full discharges (or discharging below 20%) except for this instance.

15.5 Servicing


Contact your authorized distributor. Please provide the serial number of the console which is located (1) on the product labels and (2) on the SYSTEM INFORMATION screen while the unit is in operation.

Warning *Do not open the device as it contains no User serviceable parts. Opening the device will void any warranty, and may cause damage to the device and may shock and injure you.*

Prior to returning the console for servicing, please EXPORT all patient assessment data per section 10 of the User Manual. Please do not delete any data sets, as they may be used to troubleshoot the console. All patient records may be deleted from the console prior to service completion.


15.6 Internal Sensor Calibration

The Chartis Console needs to be calibrated once every 72 months to ensure accurate sensor function. Starting one month prior to the calibration due date, the console will display this reminder message every time the device is powered up. Press **OK** to dismiss this message.

The system will display the calibration due date (MM/DD/YYYY) on the System Information screen. Calibration is performed by the manufacturer. Contact your authorized distributor to schedule a calibration service.

IMPORTANT: After the calibration due date, the Chartis Console will no longer perform assessments. You must have the console calibration performed before you will be able to assess a patient.

NOTE: Users will still be able to review stored patient data after the calibration due date.

16. Specifications

This section describes the product				
<u>Specification</u>	<u>Description</u>			
Display	1920 x 1080 Resolution 13.3 inch diagonal color touchscreen			
Approximate Dimensions (Touchscreen tablet docked and open)	Height	10.8 inches / 27.4 cm		
	Width	14.0 inches / 35.6 cm		
	Depth	11.2 inches / 28.4 cm		
Approximate Weight (Touchscreen Tablet Computer only)	3.85 lb / 1.75 kg			
Approximate Weight (Patient Interface Box only)	7.80 lb / 3.54 kg			
Input specifications	Input air flow	Range	±1500 ml/min	
		Resolution	2ml/min	
		Accuracy	±5% or ±5ml/min whichever is greater	
	Input air pressure	Range	±50 cmH ₂ O	
		Resolution	0.1 cmH ₂ O	
		Accuracy	±5% or ±0.05 cmH ₂ O whichever is greater	
	Expired air flow volume during assessment	Range	0-2 L	
		Resolution	2 ml	
		Accuracy	±5% or ±5 ml whichever is greater	
		Voltage	19 VDC ± 5%	
	Current	3.43 A ± 5%		
Operational environment	Temperature	+15°C TO +35°C (+59 TO +95 °F)		
	Humidity	15 TO 65% humidity, non-condensing		
	Maximum Altitude	2000 meter		
Storage and Transportation Environment	Temperature	-18°C TO +55°C (0 TO +131 °F)		
	Humidity	15 TO 65% humidity, non-condensing		

Device Classification	Description
Applied part type	Type BF
Medical device	Class II Equipment (IIa EU)
Electrical protection	Class I Equipment
Ingress protection	IPX0
Mode of operation	Continuous

Power supply specifications	Description
Input rating	100-240 VAC, 50-60 Hz, 1.6-0.7A
Output rating	19 VDC, 3.43 A

17. Electromagnetic Compatibility (EMC) Compliance Information

This section describes the electromagnetic emission and immunity compliance information and recommended environments. The Chartis Pulmonary Assessment System does not have any performance of a clinical function, other than that related to basic safety, that if not achieved leads to an unacceptable patient or user risk. The Chartis Pulmonary Assessment System is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment. The Chartis Pulmonary Assessment System was tested in the electromagnetic environment specified below to ensure that it functioned as intended without any unrecoverable anomalies and could shut down safely.

Emission standard	Compliance	Guidance
RF Emissions CISPR 11	Group 1	The Chartis Pulmonary Assessment System uses RF (radio frequency) energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The Chartis Pulmonary Assessment System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions EN 61000-3-2	Class A	
Voltage Fluctuations and Flicker EN 61000-3-3	Complies	

Immunity Test	IEC 60601-1-2 Test level/Compliance Level	Guidance
IEC 61000-4-2: Electrostatic Discharge (ESD)	±8 kV Contact ±15 kV Air	The relative humidity should be at least 5%.
IEC 61000-4-4: Electrical Fast Transient/burst (EFT)	±2 kV Mains	Mains power quality should be that of a typical hospital environment.
IEC 61000-4-5: Surge	±1 kV Differential Mode ±2 kV Common Mode	Mains power quality should be that of a typical hospital environment.
IEC 61000-4-8: Power Frequency (50/60 Hz) Magnetic Field	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical hospital environment.
IEC 61000-4-11: Voltage Dips/Interruptions	>95% Dip for 0.5 Cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° >95% Dip for 1 Cycle at 0°C 30% Dip for 25/30 Cycles at 0°C >95% Dip for 250/300 Cycles	Mains power quality should be that of a typical hospital environment. If the user of the Chartis Pulmonary Assessment System requires continued operation during power mains interruptions, it is recommended that the Chartis Pulmonary Assessment System be powered from an uninterruptible power supply or battery.
IEC 61000-4-39: Proximity Magnetic Fields	65 A/m, PM at 2.1 kHz, 50% duty cycle, 134.2 kHz; 7.5 A/m, PM at 50 kHz, 50% duty cycle, 13.56 MHz	Proximity magnetic fields should be that of a typical hospital environment.
IEC 61000-4-6: Conducted RF	3 Vrms 0.15 MHz – 80 MHz (6 Vrms in ISM Bands between 0.15 MHz and 80 MHz) 80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Chartis Pulmonary Assessment System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.
IEC 61000-4-3: Radiated RF	3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	

CHARTIS® CONSOLE USER MANUAL

The Chartis Pulmonary Assessment System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Chartis Pulmonary Assessment System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Chartis Pulmonary Assessment System as recommended below, according to the maximum output power of the communications equipment.


Recommended separation distances between portable and mobile RF communications equipment and the Chartis Pulmonary Assessment System:			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23



For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.


Immunity to RF Wireless Communications Equipment						
Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	Immunity Test Level (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710 745 780	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
810 870 930	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
1720 1845 1970	1700 – 1900	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation 217 Hz	2	0.3	28
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240 5500 5785	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9

18. Troubleshooting

This section describes troubleshooting symptoms and possible solutions.

Issue	Solution
<p>Touchscreen Tablet Computer not charging</p>	<p><u>Charging the Tablet Computer via the Provided External Power Supply:</u></p> <p>(a) Plug the power supply directly to the touchscreen tablet computer. This should trigger a blue LED to blink slowly on the left edge of the tablet display where a battery icon is located . Also, verify the lightning bolt symbol  is displayed next to the battery symbol on the Chartis application.</p> <div data-bbox="769 611 1138 1077" data-label="Image">  </div> <p>(b) If a solid brown light is displayed on the left edge of the tablet display where a battery icon is located , the tablet battery is dead and must be replaced. Skip to step (e).</p> <p>(c) If plugging the power supply to the tablet does not trigger a blinking blue LED, the tablet battery may have been drained substantially. Unplug and re-plug the power supply to the tablet computer a few times to activate the charging process.</p> <p>(d) Monitor the battery charge level for a few hours, and verify that the battery charge level is increasing.</p> <p>(e) Contact your authorized distributor if issue is not resolved after performing the steps above.</p> <p><u>Charging the Tablet Computer via Dock on the Patient Interface Box (PIB):</u></p> <p>(a) With the power switch on the PIB pressed to the ON position (depressed inward), verify the lightning bolt symbol  is displayed next to the battery symbol on the Chartis application.</p> <p>(b) Monitor the battery charge level for a few hours, and verify that the battery charge level is increasing.</p> <p>(c) Contact your authorized distributor if above steps were not as described.</p>

Issue	Solution
<p>Touchscreen Tablet Computer not powering on, power locked-up system, or black screen</p>	<p><u>Method 1 – Turning the tablet ON for the first time out of the box:</u></p> <p><i>Connect 19V power to the tablet computer via the PIB</i></p> <p>(a) Dock the tablet computer into the docking station of the Patient Interface Box (PIB). (b) Connect the PIB to a power outlet using the power supply provided. (c) Turn ON the PIB by pressing the power button on the front panel.</p> <p style="text-align: center;">OR</p> <p><i>Connect 19V power to the tablet computer via the power supply</i></p> <p>(a) Connect the tablet computer to a power outlet using the power supply provided.</p> <p><u>Method 2 – Force a restart:</u></p> <p>(a) Disconnect the power supply from the tablet computer if connected. (b) Press and hold the tablet’s power button for 6+ seconds. (c) The tablet computer should restart. If power is not restored, try the method below.</p> <p><i>Note: This form of hardware shutdown is intended only as a means of recovery from lockups. Do not use for normal operation.</i></p> <p><u>Method 3 – Remove all power sources:</u></p> <p>(a) Disconnect the power supply from the tablet computer if connected. (b) Remove the battery from the back side of the tablet computer.</p> <div data-bbox="527 1037 1451 1247" style="border: 1px solid black; padding: 5px;">  <p style="text-align: center;">Step 1: Push the battery latch lock</p> <p style="text-align: right;">Step 2: Slide the battery latch to the left to remove the battery from the tablet</p> </div> <p>(c) After 10 seconds, install the tablet battery and connect the power supply to the tablet computer.</p> <div data-bbox="561 1335 1333 1600" style="border: 1px solid black; padding: 5px;">  <p style="text-align: center;">Slide in the bottom of the battery first into the battery compartment followed by pushing down on the top of the battery until you hear a click.</p> </div> <p>(d) If the tablet computer does not automatically power on, press and release the power button on the tablet to power it ON.</p> <p><u>Note:</u> If tablet computer is not powering on after power cycling, wait 10+ seconds after tablet computer is OFF before attempting to power ON the tablet computer via the automatic power on methods described in section 14.2.</p> <p>Contact your authorized distributor if issue is not resolved after performing the steps above.</p>

Issue	Solution
<p>Power switch (on the Patient Interface Box) is pressed to the ON position (depressed inward) but the device doesn't power up.</p>	<p>(a) Verify that the power cable is properly plugged into the wall AC outlet on one end and the power supply on the other end. The green indicator light on the power supply should be lit.</p>  <p>(b) Verify that the barrel plug of the power cable is plugged properly into the connector at the back of the Patient Interface Box.</p> <p>(c) Confirm that the power switch on the Patient Interface Box is set to the ON position (depressed inward).</p> <p>(d) Try steps (a) through (c) with an alternate wall AC outlet.</p> <p>(e) Verify that the power supply is the unit that shipped with the device or is an approved model.</p> <p>(f) Contact your authorized distributor if issue is not resolved after performing the steps above.</p>
<p>Sensor initialization failed</p>	<ol style="list-style-type: none"> 1. Verify that both dust plugs are unplugged. 2. Verify that there is no source of disturbance near the quick disconnect socket on the back panel of the Patient Interface Box.
<p>Calibration due</p>	<p>Contact your authorized distributor.</p>
<p>Unable to connect the Chartis Catheter or reattach the dust plug</p>	<p>Push the latch of the front/back quick-disconnect socket towards the orifice and try again.</p> 
<p>Unable to detect a USB storage device</p>	<ol style="list-style-type: none"> 1. Verify that the USB device is plugged in correctly. 2. Verify that the USB device is a storage device.
<p>No waveforms are visible on the assessment screen</p>	<p>Ensure that the Patient Interface Box is powered ON.</p> 

19. Unit Conversion

This section contains the pressure conversion factor from SI pressure units, pascals (Pa), to cmH₂O (4°C).

PASCALS	cmH ₂ O (4°C)
1.00	0.010197162
98.0665	1.00

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
The screenshot displays the 'Patient Review' interface. At the top, there is a header with the 'pulmonX' logo, user information 'User TEST-USER', the title 'Patient Review', and system status including a 41% battery icon, an American flag, the date '04/12/2023', and the time '5:47 PM'. Below the header is a table with the following columns: 'Patient ID', 'Assess Date', and 'Count'. The table contains two rows: one for 'TEST-PATIENT' with an assess date of '04/12/2023 3:54 PM' and a count of '10', and another for 'TEST' with an assess date of 'None' and a count of '0'. To the right of the table are three buttons: 'Assess', 'Static review' (highlighted with a yellow circle), and 'Dynamic review'. To the right of these buttons is a 'Video (1)' button. At the bottom of the interface are two buttons: 'Export' and 'Delete'.

Patient ID	Assess Date	Count	Assess	Static review	Dynamic review	Video (1)
TEST-PATIENT	04/12/2023 3:54 PM	10				
TEST	None	0				

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The screenshot displays the 'Patient Review' interface. At the top, the user is identified as 'User TEST-USER'. The page title is 'Patient Review'. The top right corner shows a battery level of 41%, a US flag, and the date and time '04/12/2023 5:47 PM'. Below the header is a table with the following data:

<input type="checkbox"/>	Patient ID	Assess Date	Count	Assess	Static review	Dynamic review	Video (1)
<input checked="" type="checkbox"/>	TEST-PATIENT	04/12/2023 3:54 PM	10	Assess	Static review	Dynamic review	Video (1)
<input type="checkbox"/>	TEST	None	0				

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The screenshot shows the 'Patient Review' screen. At the top, there is a header with the pulmonX logo, user information (TEST-USER), patient information (TEST-PATIENT), and the title 'Patient Review'. On the right, there is a status bar showing 38% battery, a US flag, and the date/time (04/12/2023 4:05 PM). Below the header is a table with the following data:

<input type="checkbox"/>	Patient ID	Assess Date	Count	Assess	Static review	Dynamic review	Video (1)
<input checked="" type="checkbox"/>	TEST-PATIENT	04/12/2023 3:54 PM	10				
<input type="checkbox"/>	TEST	None	0				

At the bottom of the screen, there are two buttons: 'Export' and 'Delete'.

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The screenshot shows the 'Video Review' screen. At the top, there is a header with the pulmonX logo, user information (TEST-USER), patient information (TEST-PATIENT), and the title 'Video Review'. On the right, there is a status bar showing 38% battery, a US flag, and the date/time (04/12/2023 4:06 PM). Below the header is a table with the following data:

<input type="checkbox"/>	File name	Recorded Date	Size (MB)	Play
<input checked="" type="checkbox"/>	TEST-PATIENT_RML-1 VEN Video.mkv	04/12/2023 4:02 PM	2.7	

At the bottom of the screen, there is a 'Delete' button.

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
The screenshot shows the 'Patient Review' interface. At the top, there is a header with the 'pulmonX' logo, 'User TEST-USER', and the title 'Patient Review'. On the right, there is a status bar showing 41% battery, a US flag, the date '04/12/2023', and the time '5:47 PM'. Below the header is a table with the following columns: 'Patient ID', 'Assess Date', and 'Count'. The table contains two rows: one for 'TEST-PATIENT' with an assess date of '04/12/2023 3:54 PM' and a count of '10', and another for 'TEST' with an assess date of 'None' and a count of '0'. To the right of the table are four buttons: 'Assess', 'Static review' (highlighted with a yellow circle), 'Dynamic review', and 'Video (1)'. At the bottom of the interface are two buttons: 'Export' and 'Delete'.

Patient ID	Assess Date	Count	Assess	Static review	Dynamic review	Video (1)
TEST-PATIENT	04/12/2023 3:54 PM	10				
TEST	None	0				

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Copyright © 2023 by Pulmonx Corporation. Alle Rechte vorbehalten.
Copyright © 2023 di Pulmonx Corporation. Tutti i diritti riservati.
Copyright © 2023 Pulmonx Corporation. Reservados todos los derechos.
Copyright © 2023 par Pulmonx Corporation. Tous droits réservés.
Copyright © 2023 by Pulmonx Corporation. Todos os direitos reservados.
Copyright © 2023 af Pulmonx Corporation. Alle rettigheder forbeholdes.